

**Model No. DPN-526-XFYP2Q**  
**526.5nm DPSS LASER UP TO >1000mW**


The DPN-series 526.5nm LD pumped all-solid-state laser is constructed with features of high stability, FDA compliance, ultra-compactness, long lifetime, and easy operating. It is widely used in collimation, laser medical treatment, scientific experiment, optical instrument, and many other applications.

**SPECIFICATIONS**

<b>Model No.</b>	<b>DPN-526-XFYP2Q</b>
<b>Wavelength (nm)</b>	526.5±1
<b>Output power (mW)</b>	>300 ( <b>X=300</b> ), >500 ( <b>X=500</b> ), >800 ( <b>X=800</b> ), >1000 ( <b>X=1W</b> )
<b>Transverse mode</b>	Near TEM <sub>00</sub>
<b>Operating mode</b>	CW
<b>Power stability (rms, over 4 hours)</b>	<5% ( <b>P=A</b> ), <3% ( <b>P=E</b> )
<b>Dimensions of beam at aperture (1/e<sup>2</sup>, mm)</b>	~3.0
<b>Beam divergence, full angle (mrad)</b>	<1.5
<b>M<sup>2</sup> factor</b>	3~5
<b>Polarization ratio</b>	>100:1
<b>Warm-up time (minutes)</b>	<10
<b>Beam height from base plate (mm)</b>	68.5
<b>Operating temperature (°C)</b>	10~35
<b>Laser head</b>	245 (L) ×99 (W) ×94 (H) mm <sup>3</sup> , 2.6 kg
<b>Power supply (90-240VAC)</b>	<b>High Power Elite Power Supply (Y=H):</b> 307 (L) ×150 (W) ×106 (H) mm <sup>3</sup> , 2.9 kg; with complete FDA compliant features, such as the turnkey switch and interlock, and easy to operate; CW mode with optional TTL modulation up to 10kHz  <b>High Power Laboratory Power Supply (Y=M):</b> 307 (L) ×150 (W) ×106 (H) mm <sup>3</sup> , 2.9 kg; with complete FDA compliant features, such as the turnkey switch and interlock, and more functions; CW mode with the adjustable output power knob, the operating current LED display, and optional TTL/Analog modulations up to 10kHz
<b>Optional Modulation</b>	None ( <b>Q=0</b> ) <b>TTL:</b> 1Hz -1kHz ( <b>Q=1</b> ), 1kHz -10kHz ( <b>Q=2</b> ) <b>Analog:</b> 1Hz -1kHz ( <b>Q=3</b> ), 1kHz -10kHz ( <b>Q=4</b> )
<b>Expected lifetime (hours)</b>	10,000
<b>Warranty</b>	10 months
<b>Remarks</b>	The stability of output power is promised at the maximum output power only. The stability of output power may change when you adjust the output power.


**Note: The above specifications are subject to change without notice.**






**Mechanical Outline (unit: mm)**

**Laser Head**




**Power Supply**

**High Power Elite Power Supply (Y=H)**



**High Power Laboratory Power Supply (Y=M)**



**Note: The above specifications are subject to change without notice.**

